

PCNA Rabbit mAb

Catalog No: #52903

Package Size: #52903-1 50ul #52903-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PCNA Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human;Mouse;Rat
Specificity	PCNA Antibody detects endogenous levels of total PCNA
Immunogen Description	A synthesized peptide derived from human PCNA
Conjugates	Unconjugated
Other Names	Cyclin; DNA polymerase delta auxiliary protein; MGC8367; PCNA; Proliferating cell nuclear antigen
Accession No.	Swiss Prot:P12004
Calculated MW	29kDa
Formulation	Liquid
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB 1:3000~1:10000

IHC 1:50~1:200 ICC/IF 1:100~1:500

IP 1:50~1:100 FC 1:200~1:500

Images



Western blot analysis of extracts from 3T3, 293, HeLa, Jurkat and MCF7 using PCNA Rabbit mAb #52903.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.