

Lamin B1 Rabbit mAb

Catalog No: #52904

Package Size: #52904-1 50ul #52904-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Lamin B1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human;Mouse;Rat
Immunogen Description	A synthesized peptide derived from human Lamin B1
Conjugates	Unconjugated
Other Names	LMN; ADLD; LMN2; LMNB; Lamin B1;
Accession No.	Swiss Prot:P20700
Calculated MW	70kDa
Formulation	Liquid
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:3000~1:10000

IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of extracts from 3T3, 293, HeLa, Mouse Brain and Rat Brain using Lamin B1 Rabbit mAb #52904.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.