

β-Tubulin Rabbit mAb

Catalog No: #52905

Package Size: #52905-1 50ul #52905-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	β-Tubulin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat Cow Dog Chicken Monkey Zebrafish Xenopus tropicalis Drosophila melanogaster
Specificity	beta Tubulin Antibody detects endogenous levels of total beta Tubulin
Immunogen Description	A synthesized peptide derived from human beta Tubulin
Other Names	Beta 4 tubulin; Beta 5 tubulin; Beta1 tubulin; Tubulin beta chain; Tubulin beta class I;
Accession No.	Swiss Prot:P07437
Calculated MW	50kDa
Formulation	Liquid
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:3000~1:10000

IHC 1:100~1:200 ICC/IF 1:100~1:500 FC 1:50

Images



Western blot analysis of extracts from 3T3, 293, HeLa, Mouse Brain and Rat Brain using beta Tubulin Rabbit mAb Antibody #52905.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.