

APOF Rabbit Polyclonal Antibody

Catalog No: #54279



Package Size: #54279-1 50ul #54279-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	APOF Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant fusion protein of human APOF (NP_001629.1).
Conjugates	Unconjugated
Other Names	APOF;Apo-F;LTIP
Accession No.	Swiss Prot:Q13790GeneID:319
Calculated MW	35kDa
Formulation	Buffer: PBS with 0.02% sodium azide, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:1000 IHC □ 1:50 - 1:100 IF □ 1:50 - 1:100

Background

The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.