

## SERINC1 Rabbit Polyclonal Antibody

Catalog No: #54550



Package Size: #54550-1 50ul #54550-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SERINC1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human;Mouse;Rat
Immunogen Description	Recombinant protein of human SERINC1.
Conjugates	Unconjugated
Other Names	TDE1L;TDE2;TMS-2;TMS2;SERINC1
Accession No.	Swiss Prot:Q9NRX5GeneID:57515
Calculated MW	–
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200

## Background

–

## Published Papers

el at., Genomic landscape defines peritoneal metastatic pattern and related target of peritoneal metastasis in colorectal cancer. In Int J Cancer on 2024 Oct 1 by Minhui Hu, Rui Luo,et al..PMID:38738976, , (2024)

[PMID:38738976](https://pubmed.ncbi.nlm.nih.gov/38738976/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.