DHFR Rabbit mAb

Catalog No: #56208

Package Size: #56208-1 50ul #56208-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com


Description

Product Name	DHFR Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human;Mouse;Rat
Specificity	DHFR Antibody detects endogenous levels of total DHFR
Immunogen Description	A synthesized peptide derived from human DHFR
Conjugates	Unconjugated
Other Names	DHFR; DHFRP1; Dihydrofolate reductase; DYR;
Accession No.	Uniprot:P00374
Calculated MW	21kDa
SDS-PAGE MW	21kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:5000 IHC:1:50~1:200 ICC/IF:1:50~1:200

Images



Western blot analysis of DHFR expression in (1) HeLa cell lysate; (2) Mouse spleen lysate.

Product Description

Key enzyme in folate metabolism. Catalyzes an essential reaction for de novo glycine and purine synthesis, and for DNA precursor synthesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.		