## NCX1 Rabbit mAb

Catalog No: #56340

Package Size: #56340-1 50ul #56340-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com


# Description

Product Name	NCX1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	NCX1 Antibody detects endogenous levels of total NCX1
Immunogen Description	A synthesized peptide derived from human NCX1
Conjugates	Unconjugated
Other Names	Na+/Ca2+ exchanger; NCX1; SLC8A1;
Accession No.	Uniprot:P32418
Calculated MW	108kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## **Application Details**

WB:1:500~1:2000

## **Images**



Western blot analysis of NCX1 expression in K562 cell lysate.

## **Product Description**

Rapidly transports Ca(2+) during excitation-contraction coupling. Ca(2+) is extruded from the cell during relaxation so as to prevent overloading of intracellular stores.

## **Published Papers**

Menglong Wang; Tongrui Wang; Yang Liu; Lurong Zhou; Yuanping Yin; Feng Gu el at., Identification and study of mood-related biomarkers and potential molecular mechanisms in type 2 diabetes mellitus., , (2025)

PMID:39915429

Note: This product is for in vitro research use only and is not intended for use in humans or animals.