GOLPH4 Rabbit mAb

Catalog No: #56511

Package Size: #56511-1 50ul #56511-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description


Product Name	GOLPH4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	GOLPH4 Antibody detects endogenous levels of total GOLPH4
Immunogen Description	A synthesized peptide derived from human GOLPH4
Conjugates	Unconjugated
Other Names	GIMPc; golim4; GOLPH 4; GPP130; P 138; P138;
Accession No.	Uniprot:O00461
Calculated MW	82kDa
SDS-PAGE MW	130kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:100~1:500

Images



Western blot analysis of GOLPH4 expression in HepG2 cell lysate.

Product Description

Plays a role in endosome to Golgi protein trafficking; mediates protein transport along the late endosome-bypass pathway from the early endosome to the Golgi.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.			