

GOLPH4 Rabbit mAb

Catalog No: #56511

Package Size: #56511-1 50ul #56511-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GOLPH4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	GOLPH4 Antibody detects endogenous levels of total GOLPH4
Immunogen Description	A synthesized peptide derived from human GOLPH4
Conjugates	Unconjugated
Other Names	GIMPC; golim4; GOLPH 4; GPP130; P 138; P138;
Accession No.	Uniprot:O00461
Calculated MW	82kDa
SDS-PAGE MW	130kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

IHC:1:100~1:500

Images



All lanes: GOLPH4 Rabbit mAb at 1/1k dilution

Lane 1 : HepG2 whole cell lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 82 kDa

Observed band size: 130 kDa

Exposure time: 6 seconds




All lanes: GOLPH4 Rabbit mAb at 1/1k dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : GOLPH4 knockdown HeLa cell lysate

Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human colon tissue stained for GOLPH4 using 56511 at 1/100 dilution in immunohistochemical analysis.

Product Description

Plays a role in endosome to Golgi protein trafficking; mediates protein transport along the late endosome-bypass pathway from the early endosome to the Golgi.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.