

CBX2 Rabbit mAb

Catalog No: #56608

Package Size: #56608-1 50ul #56608-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CBX2 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	CBX2 Antibody detects endogenous levels of total CBX2
Immunogen Description	A synthesized peptide derived from human CBX2
Conjugates	Unconjugated
Other Names	CBX2; CDCA6; M33 polycomb like; MOD2; Modifier3; Polycomb; SRXY5;
Accession No.	Q14781
Calculated MW	56 kDa
SDS-PAGE MW	56 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:1000~1:5000

Images



All lanes : CBX2 Rabbit mAb at 1/1k dilution Lane 1 : HeLa whole cell lysates Lane 2 : SH-SY5Y whole cell lysates Lane 3 : MCF7 whole cell lysates Lane 4 : 3T3 whole cell lysates Lane 5 : C6 whole cell lysates Lane 6 : Rat Brain lysates Lysates/proteins at 20 µg per lane. Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution Predicted band size: 56 kDa Observed band size: 56 kDa Exposure time: 1 seconds



All lanes : CBX2 Rabbit mAb at 1/1k dilution Lane 1 : Wild-type HeLa cell lysate Lane 2 : CBX2 knockdown HeLa cell lysate Lysates/proteins at 20 µg per lane.

Product Description

Component of the Polycomb group (PcG) multiprotein PRC1 complex, a complex required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.