

HEAB Rabbit mAb

Catalog No: #56992



Package Size: #56992-1 50ul #56992-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HEAB Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	HEAB Antibody detects endogenous levels of total HEAB
Immunogen Description	A synthesized peptide derived from human HEAB
Conjugates	Unconjugated
Other Names	clp1; hClp1;
Accession No.	Uniprot:Q92989
Calculated MW	48kDa
SDS-PAGE MW	48kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000

IHC:1:50~1:200

FC:1:50

Images



All lanes: HEAB Rabbit mAb at 1/1k dilution

Lane 1 : Hela whole cell lysates
Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 48 kDa

Observed band size: 48 kDa

Exposure time: 10 seconds




All lanes:HEAB Rabbit mAb at 1/1k dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : HEAB knockdown HeLa cell lysate

Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human colon cancer tissue stained for HEAB using 56992 at 1/100 dilution in immunohistochemical analysis.

Product Description

Polynucleotide kinase that can phosphorylate the 5'-hydroxyl groups of double-stranded RNA (dsRNA), single-stranded RNA (ssRNA), double stranded DNA (dsDNA) and double-stranded DNA:RNA hybrids.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.