

SEC62 Rabbit mAb

Catalog No: #57144



Package Size: #57144-1 50ul #57144-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SEC62 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human;Mouse;Rat
Specificity	SEC62 Antibody detects endogenous levels of total SEC62
Immunogen Description	A synthesized peptide derived from human SEC62
Conjugates	Unconjugated
Other Names	Mouse, Rat, Human
Accession No.	Q99442
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of SEC62 expression in HeLa cell lysate.

Background

Required for preprotein translocation.

Published Papers

Zhe Chen;Jun Cheng;Qun Zhou;Le-le Wu;Jia-Wei Chen;Xiang-Ning Duan;Jia-Long Yan;Jian-Gang Cao;Xiao-Dan Xia;Lan-Fang Li;Lin-Xi Chen et al.,
SEC62-dependent ER-phagy contributes to apelin-13/APJ-induced monocyte-vascular endothelial cell adhesion in atherosclerosis pathogenesis., ,
(2025)

[PMID:39930135](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.