

Cyclooxygenase 2 Rabbit mAb

Catalog No: #57246



Package Size: #57246-1 50ul #57246-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cyclooxygenase 2 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of mouse COX2
Conjugates	Unconjugated
Other Names	Cyclooxygenase-2; COX-2; Glucocorticoid-regulated inflammatory cyclooxygenase; Grippghs; PES-2; PHS II; Prostaglandin H2 synthase 2
Accession No.	Swiss-Prot:Q05769 Gene ID:19225
Calculated MW	69 kDa
SDS-PAGE MW	74 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

Images



Western blot analysis of COX2 in 3t3 lysates using Cyclooxygenase 2 antibody.

Published Papers

Xia Zhang et al., analysis: an analysis based on the European spontaneous adverse event reporting system, *Frontiers in Pharmacology*, (2025)

PMID:39877384

Note: This product is for in vitro research use only and is not intended for use in humans or animals.