

Hsp90 beta Rabbit mAb

Catalog No: #58631



Package Size: #58631-1 50ul #58631-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hsp90 beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Hsp90 beta Antibody detects endogenous levels of total Hsp90 beta
Conjugates	Unconjugated
Other Names	HS90B;HSP 84;HSP90-beta;HSP90AB1; HSPC2; HSPCB; Heat shock protein HSP 90-beta;
Accession No.	P08238
Calculated MW	83 kDa
SDS-PAGE MW	90 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000


IHC: 1:50-1:200

ICC/IF: 1:50-1:200


Images




All lanes: Hsp90 beta Rabbit mAb at 1/1k dilution
 Lane 1 : HeLa whole cell lysates
 Lane 2 : JK whole cell lysates
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 83 kDa
 Observed band size: 90 kDa
 Exposure time: 7 seconds




All lanes: Hsp90 beta Rabbit mAb at 1/1k dilution
 Lane 1 : Mouse heart lysates
 Lane 2 : Mouse liver lysates
 Lane 3 : Mouse brain lysates
 Lane 4 : Rat heart lysates
 Lane 5 : Rat liver lysates
 Lane 6 : Rat brain lysates
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 83 kDa
 Observed band size: 90 kDa
 Exposure time: 6 seconds



All lanes :Hsp90 beta Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HAP1 cell lysate
 Lane 2 : Hsp90 beta knockdown HAP1 cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human testis tissue stained for Hsp90 beta using 58631 at 1/100 dilution in immunohistochemical analysis.



Immunocytochemistry/ Immunofluorescence Hsp90 beta antibody (58631) ICC/IF staining of Hsp90 beta in HeLa cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 58631 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 488 goat anti rabbit, used at a dilution of 1/500.

Nuclei were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.