

JAM1 Rabbit mAb

Catalog No: #58918



Package Size: #58918-1 50ul #58918-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	JAM1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	JAM1 Antibody detects endogenous levels of JAM1
Conjugates	Unconjugated
Other Names	F11R;CD321;JAM;JAM-1;JAM-A;JAM1;JAMA;JCAM;KAT;PAM-1;Junction adhesion molecule 1;
Accession No.	Q9Y624
Calculated MW	33 kDa
SDS-PAGE MW	33 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000


IHC: 1:50-1:200

ICC/IF: 1:50-1:200


Images




All lanes: JAM1 Rabbit mAb at 1/1k dilution
 Lane 1 : HeLa whole cell lysates Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 33 kDa Observed band size: 33 kDa
 Exposure time: 4 seconds



All lanes: JAM1 Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HeLa cell lysate
 Lane 2 : JAM1 knockdown HeLa cell lysate
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human bladder cancer tissue stained for JAM1 using 58918 at 1/100 dilution in immunohistochemical analysis.



Immunocytochemistry/ Immunofluorescence JAM1 antibody (58918)
 ICC/IF staining of JAM1 in MCF7 cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 58918 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 488 goat anti rabbit, used at a dilution of 1/500.

Nuclei were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.