

## PRDX1 Rabbit mAb

Catalog No: #59228

Package Size: #59228-1 50ul #59228-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	PRDX1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;IHC;ICC/IF
Species Reactivity	Human Mouse
Specificity	PRDX1 Antibody detects endogenous levels of total PRDX1
Conjugates	Unconjugated
Other Names	PAG; PAGA; PAGB; PRX1; PRXI; MSP23; NKEFA; TDPX2; NKEF-A; Peroxiredoxin1; PRDX1;
Accession No.	Q06830
Calculated MW	22 kDa
SDS-PAGE MW	22 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


## Application Details

WB: 1:500-1:2000


IHC: 1:50-1:200

ICC/IF: 1:50-1:200


## Images




All lanes: PRDX1 Rabbit mAb at 1/1k dilution  
 Lane 1 : HeLa whole cell lysates  
 Lane 2 : JK whole cell lysates  
 Lane 3 : HepG2 whole cell lysates  
 Lane 4 : MCF-7 whole cell lysates  
 Lysates/proteins at 20 µg per lane.  
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution  
 Predicted band size: 22 kDa  
 Observed band size: 22 kDa  
 Exposure time: 4 seconds




All lanes: PRDX1 Rabbit mAb at 1/1k dilution  
 Lane 1 : Mouse kidney lysates  
 Lane 2 : Mouse brain lysates  
 Lane 3 : Mouse skin lysates  
 Lysates/proteins at 20 µg per lane.  
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution  
 Predicted band size: 22 kDa  
 Observed band size: 22 kDa  
 Exposure time: 5 seconds



All lanes: PRDX1 Rabbit mAb at 1/1k dilution  
 Lane 1 : Wild-type HAP1 cell lysate  
 Lane 2 : PRDX1 Rabbit mAb knockdown HAP1 cell lysate  
 Lysates/proteins at 20 µg per lane.



Formalin-fixed, paraffin-embedded human thyroid tissue stained for PRDX1 using 59228 at 1/100 dilution in immunohistochemical analysis.



Immunocytochemistry/ Immunofluorescence PRDX1 antibody (59228)  
 ICC/IF staining of PRDX1 in 293T cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 59228 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 647 goat anti rabbit, used at a dilution of 1/500.

The negative control is shown in bottom right hand panel - for the negative control. Nuclei were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.