

ATG7 Rabbit mAb

Catalog No: #59549



Package Size: #59549-1 50ul #59549-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ATG7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;ICC/IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	A synthesized peptide derived from human ATG7
Conjugates	Unconjugated
Other Names	hAGP7; Ubiquitin-activating enzyme E1-like protein; APG7L
Accession No.	O95352
Calculated MW	78 kDa
SDS-PAGE MW	77 kDa
Formulation	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.


Application Details

WB 1:1000-1:2000 ICC/IF 1:50-1:200


Images




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes : ATG7 Rabbit mAb at 1/1k dilution
Lane 1 : Wild-type HAP1 cell lysate
Lane 2 : ATG7 knockdown HAP1 cell lysate
Lysates/proteins at 20 µg per lane.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.