

AAMP Rabbit mAb

Catalog No: #60181



Package Size: #60181-1 50ul #60181-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	AAMP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB;ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	AAMP Antibody detects endogenous levels of total AAMP
Conjugates	Unconjugated
Other Names	AAMP;
Accession No.	Q13685
Calculated MW	47 kDa
SDS-PAGE MW	52 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

ICC/IF: 1:50-1:200


Images



All lanes: AAMP Rabbit mAb at 1/1k dilution
Lane 1 : A375 whole cell lysates
Lane 2 : MCF-7 whole cell lysates
Lane 3 : Mouse brain lysates
Lane 4 : Rat brain lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 47 kDa
Observed band size: 52 kDa
Exposure time: 6 seconds



All lanes: AAMP Rabbit mAb at 1/1k dilution
 Lane 1 : Wild-type HeLa cell lysate
 Lane 2 : AAMP Rabbit mAb knockdown HeLa cell lysate
 Lysates/proteins at 20 µg per lane.



Immunocytochemistry/ Immunofluorescence AAMP antibody (60181)
 ICC/IF staining of AAMP in HeLa cells. Cells were fixed with 4% Paraformaldehyde permeabilized with 0.1% Triton X-100.

Samples were incubated with 60181 at a working dilution of 1/100. The secondary antibody was Alexa FluorB 647 goat anti rabbit, used at a dilution of 1/500.

The negative control is shown in bottom right hand panel - for the negative control. Nuclei were counterstained with DAPI.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.