

Rat anti-Mouse CD3 Monoclonal Antibody

Catalog No: #64440

Package Size: #64440-1 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat anti-Mouse CD3 Monoclonal Antibody
Host Species	Rat
Clonality	Monoclonal
Isotype	Rat IgG2a
Purification	protein A purified
Applications	FCM
Species Reactivity	Mouse
Conjugates	Unconjugated
Target Name	Cd3e
Other Names	T cell antigen receptor complex, T3
Accession No.	P22646
Formulation	Purified antibody in PBS with 0.05% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

FCM 1:100-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.