

Wnt3 Mouse mAb

Catalog No: #64441

Package Size: #64441-1 50ul #64441-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Wnt3 Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Transfected
Specificity	Transfected Only.
Immunogen Description	Purified recombinant human Wnt3 protein fragments expressed in E.coli.
Conjugates	Unconjugated
Target Name	wnt3
Other Names	Int 4 ; INT4 ; MGC131950 ; MGC138321 ; MGC138323 ; Proto-oncogene Int-4 homolog ; Proto-oncogene Wnt-3 ; wingless type MMTV integration site family member 3 ; Wnt 3 proto oncogene protein ; WNT 3 proto oncogene protein precursor ; wnt3 ; WNT3_HUMAN.
Accession No.	P56703
Calculated MW	40kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000

Images



Western blot detection of Wnt3 in CHO-K1 cell lysate (B) and CHO-K1 transfected by Wnt3 (A) cell lysate using Wnt3 mouse mAb (1:1000 diluted). Predicted band size: 40KDa. Observed band size: 45KDa.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.