

Recombinant Human High Mobility Group Protein B3, HMGB3 (C-6His)

Catalog No: #AP75401



Package Size: #AP75401-1 10ug #AP75401-2 50ug #AP75401-3 500ug #AP75401-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human High Mobility Group Protein B3, HMGB3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O15347
Calculated MW	24kDa
Target Sequence	MAKGDPKKPKGKMSAYAFFVQTCREEHKKKNPEVPVNFAEFSKKCSERWKTMSGKEKSKFDEMAKADKVR YDREMKDYGPAKGGKKKKDPNAPKRPPSGFFLCSEFRPKIKSTNPGISIGDVAKKLGEMWNNLNDSEKQPYI TKAAKLKEYEKDVADYKSKGKFDGAKGPAKVARKKVEEEDEEEEEEEEEEEEEDEVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.