

Recombinant Human Transforming Growth Factor β -1,TGFB1

Catalog No: #AP75479

Package Size: #AP75479-1 10ug #AP75479-2 50ug #AP75479-3 500ug #AP75479-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Transforming Growth Factor β -1,TGFB1
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01137
Calculated MW	12.8kDa
Target Sequence	ALDTNYCFSSSTEKNCCVRQLYIDFRKDLGWKWIHEPKGYHANFCLGPCPYIWSLDTQYSKVLALYNQHNP GA SAAPCCVPQALEPLPIVYYYVGRKPKVEQLSNMIVRSCKCS
Formulation	Lyophilized from a 0.2 μ m filtered solution of 50mM Glycine 50mM NaCl pH4.0.
Storage	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7°C for 2-7 days.Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Published Papers

el at., Comprehensive Prognostic Analysis of Immune Implication Value and Oxidative Stress Significance of NECAP2 in Low-Grade Glioma.In Oxid Med Cell Longev. 2022 Dec 7 by Zhichao Lu, Yixun Chen , et al..PMID:36531205, , (2022)

[PMID:36531205](https://pubmed.ncbi.nlm.nih.gov/36531205/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.