

Recombinant Mouse GM-CSF,CSF2

Catalog No: #AP76135

Package Size: #AP76135-1 10ug #AP76135-2 50ug #AP76135-3 500ug #AP76135-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse GM-CSF,CSF2
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P01587
Calculated MW	14.2kDa
Target Sequence	MAPTRSPITVTRPWKHVEAIKEALNLLDDMPVTLNEEVEVVSNEFSFKKLTVCVQTRLKIFEQGLRGNFTKLKGA LNMTASYQTYCPPTPETDCETQVTYADFIDSLKTFLTDPFECKKPGQK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7°C for 2-7 days.Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Published Papers

el at., Multiple mechanisms involved in diabetes protection by lipopolysaccharide in non-obese diabetic mice. In Toxicol Appl Pharmacol
on 2015 Jun 15 by Jun Wang, Hui Cao,et al..PMID:25896969, , (2015)

[PMID:25896969](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.