

Recombinant Human FGF2 Protein, N-His-SUMO

Catalog No: #AP90760

Package Size: #AP90760-1 50ug #AP90760-2 100ug #AP90760-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FGF2 Protein, N-His-SUMO
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	bFGF, FGF2, Basic fibroblast growth factor, HBGF-2, FGFB, FGF-2, Heparin-binding growth factor 2, Fibroblast growth factor 2
Accession No.	P09038
Target Species	Homo sapiens (Human)
Calculated MW	28.66 kDa
Target Length	Pro143-Ser288
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.