

## Recombinant Human RSPO1 Protein, N-His

Catalog No: #AP92408



Package Size: #AP92408-1 50ug #AP92408-2 100ug #AP92408-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RSPO1 Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	R-spondin-1, hRspo1, RSPO1, Roof plate-specific spondin-1
Accession No.	Q2MKA7
Target Species	Homo sapiens (Human)
Calculated MW	15.00 kDa
Target Length	Ser21-Ala135
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.