

Recombinant Mouse NOG/Noggin Protein, C-Fc

Catalog No: #AP95254



Package Size: #AP95254-1 50ug #AP95254-2 100ug #AP95254-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse NOG/Noggin Protein, C-Fc
Host Species	Mammalian Cells
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Noggin, Nog
Accession No.	P97466
Target Species	Mus musculus (Mouse)
Calculated MW	53 kDa
Target Length	Met1-Cys232
Tag Info	Fc-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.