

Recombinant Mouse FGF2/bFGF Protein, N-His

Catalog No: #AP95826

Package Size: #AP95826-1 50ug #AP95826-2 100ug #AP95826-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse FGF2/bFGF Protein, N-His
Host Species	E. coli
Purification	>90% as determined by SDS-PAGE.
Applications	ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress
Other Names	Fibroblast growth factor 2, Basic fibroblast growth factor, Heparin-binding growth factor 2, Fgf2, Fgf-2, HBGF-2, FGF-2, bFGF
Accession No.	P15655
Target Species	Mus musculus (Mouse)
Calculated MW	18.50 kDa
Target Length	Ala11-Ser154
Tag Info	His-tag
Formulation	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.