

Recombinant Mouse EGF Protein, C-Fc

Catalog No: #AP96077



Package Size: #AP96077-1 50ug #AP96077-2 100ug #AP96077-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|----------------|---|
| Product Name | Recombinant Mouse EGF Protein, C-Fc |
| Host Species | Mammalian Cells |
| Purification | >90% as determined by SDS-PAGE. |
| Applications | ELISA, Immunogen, SDS-PAGE, WB, Bioactivity testing in progress |
| Other Names | Pro-epidermal growth factor, EGF, Epidermal growth factor, Egf |
| Accession No. | P01132 |
| Target Species | Mus musculus (Mouse) |
| Calculated MW | 37.12 kDa |
| Target Length | Asn977-Arg1029 |
| Tag Info | Fc-tag |
| Formulation | Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol. |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.