

alpha Synuclein Antibody FITC Conjugated

Catalog No: #C00005F

Package Size: #C00005F 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	alpha Synuclein Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Ra
Immunogen Description	KLH conjugated synthetic peptide aa 45-95 140 derived from human alpha Synuclein
Conjugates	FITC
Target Name	alpha Synuclein
Other Names	PD1; NACP; PARK1; PARK4; Alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; SNCA
Cell Localization	Intracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

May be involved in the regulation of dopamine release and transport. Induces fibrillization of microtubule-associated protein tau. Reduces neuronal responsiveness to various apoptotic stimuli, leading to a decreased caspase-3 activation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.