

Gastrin Releasing Peptide Antibody HRP Conjugated

Catalog No: #C00006H

Package Size: #C00006H 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gastrin Releasing Peptide Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide derived from human GRP
Conjugates	HRP
Target Name	Gastrin Releasing Peptide
Other Names	BN; GRP-10; proGRP; preproGRP; Gastrin-releasing peptide; GRP
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

GRP stimulates gastrin release as well as other gastrointestinal hormones. Operates as a negative feedback regulating fear and established a causal relationship between GRP-receptor gene expression, long-term potentiation, and amygdala-dependent memory for fear (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.