

Osteopontin Antibody Biotin Conjugated

Catalog No: #C00012B

Package Size: #C00012B 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Osteopontin Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 150-200 314 derived from human Osteopontin
Conjugates	Biotin
Target Name	Osteopontin
Other Names	Osteopontin; Bone sialoprotein 1; Nephropontin; Secreted phosphoprotein 1; SPP-1; Urinary stone protein; Uropontin; SPP1; BNSP; OPN; PSEC0156
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

Background

Acts as a cytokine involved in enhancing production of interferon-gamma and interleukin-12 and reducing production of interleukin-10 and is essential in the pathway that leads to type I immunity (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.