

TGF alpha Antibody PE Conjugated

Catalog No: #C00049P

Package Size: #C00049P 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGF alpha Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	FCM
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 24-52 160 derived from human TGF-alpha
Conjugates	PE
Target Name	TGF alpha
Other Names	TFGA; Protransforming growth factor alpha; TGFA
Accession No.	Swiss-Prot#P01135NCBI Gene ID7039
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

FCM: 1:20-100

Background

TGF alpha is a mitogenic polypeptide that is able to bind to the EGF receptor EGFR and to act synergistically with TGF beta to promote anchorage-independent cell proliferation in soft agar.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.