

# NKB Antibody Biotin Conjugated

Catalog No: #C00052B

Package Size: #C00052B 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NKB Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 81-94 121 derived from human NKB
Conjugates	Biotin
Target Name	NKB
Other Names	NKB; HH1; NKNB; PRO1155; ZNEUROK1; Tachykinin-3; TAC3; UNQ585 PRO1155
Accession No.	Swiss-Prot#Q9UHF0NCBI Gene ID6866
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

Immunohistochemistry1:100-500

## Background

Tachykinins are active peptides which excite neurons, evoke behavioral responses, are potent vasodilators and secretagogues, and contract (directly or indirectly) many smooth muscles (By similarity). Is a critical central regulator of gonadal function.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.