

Neuropeptide Y Antibody Biotin Conjugated

Catalog No: #C00053B

Package Size: #C00053B 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuropeptide Y Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 36-48 97 derived from human NPY
Conjugates	Biotin
Target Name	Neuropeptide Y
Other Names	PYY4; Pro-neuropeptide Y; NPY
Accession No.	Swiss-Prot#P01303NCBI Gene ID4852
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Immunohistochemistry1:100-500

Background

NPY is implicated in the control of feeding and in secretion of gonadotrophin-release hormone.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.