

# Neuropeptide Y Antibody HRP Conjugated

Catalog No: #C00053H

Package Size: #C00053H 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                 |
|-----------------------|---------------------------------------------------------------------------------|
| Product Name          | Neuropeptide Y Antibody HRP Conjugated                                          |
| Host Species          | Rabbit                                                                          |
| Clonality             | Polyclonal                                                                      |
| Isotype               | IgG                                                                             |
| Purification          | Purified by Protein A.                                                          |
| Applications          | IHC-P                                                                           |
| Species Reactivity    | Hu Ms Rt                                                                        |
| Immunogen Description | KLH conjugated synthetic peptide aa 36-48 97 derived from human NPY             |
| Conjugates            | HRP                                                                             |
| Target Name           | Neuropeptide Y                                                                  |
| Other Names           | PYY4; Pro-neuropeptide Y; NPY                                                   |
| Accession No.         | Swiss-Prot#P01303NCBI Gene ID4852                                               |
| Cell Localization     | Secreted                                                                        |
| Concentration         | 1mg ml                                                                          |
| Formulation           | 10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% Gentamicin. |
| Storage               | Store at 4C for 12 months.                                                      |

## Application Details

Immunohistochemistry1:100-500

## Background

NPY is implicated in the control of feeding and in secretion of gonadotrophin-release hormone.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.