

Neurturin Antibody Biotin Conjugated

Catalog No: #C00055B

Package Size: #C00055B 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurturin Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 140-197 197 derived from human Neurturin
Conjugates	Biotin
Target Name	Neurturin
Other Names	NTN; Neurturin; NRTN
Accession No.	Swiss-Prot#Q99748NCBI Gene ID4902
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

Background

Supports the survival of sympathetic neurons in culture. May regulate the development and maintenance of the CNS. Might control the size of non-neuronal cell population such as haemopoietic cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.