

Nociceptin Antibody Biotin Conjugated

Catalog No: #C00057B

Package Size: #C00057B 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nociceptin Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 130-142 176 derived from human Nociceptin
Conjugates	Biotin
Target Name	Nociceptin
Other Names	PPNOC; Prepronociceptin; PNOC; OFQ
Accession No.	Swiss-Prot#Q13519NCBI Gene ID5368
Cell Localization	Secreted
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

Western blotting: 1:100-1000Immunohistochemistry1:100-500

Background

Nociceptin is the ligand of the opioid receptor-like receptor (OPRL1). It may act as a transmitter in the brain by modulating nociceptive and locomotor behavior. May be involved in neuronal differentiation and development (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.