

# Secretin receptor Antibody FITC Conjugated

Catalog No: #C00067F

Package Size: #C00067F 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Secretin receptor Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 170-220 440 derived from human of human Secretin receptor
Conjugates	FITC
Target Name	Secretin receptor
Other Names	SR; Secretin receptor; SCT-R; SCTR
Accession No.	Swiss-Prot#P47872NCBI Gene ID6344
Cell Localization	Helical, Transmembrane
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

Immunofluorescence: 1:50-200

## Background

This is a receptor for secretin. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.