

Recombinant Rabbit IgG Isotype Control Conjugated Antibody

Catalog No: #CIC004



Orders: order@signalwayantibody.com


Support: tech@signalwayantibody.com

Package Size: #CIC004-AF350 100ul #CIC004-AF405 100ul #CIC004-AF488 100ul #CIC004-AF555 100ul #CIC004-AF594 100ul
#CIC004-AF647 100ul #CIC004-AF680 100ul #CIC004-AF750 100ul #CIC004-Biotin 100ul #CIC004-Conjugated 50ul

Description

Product Name	Recombinant Rabbit IgG Isotype Control Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Formulation	Liquid in PBS containing 40% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months.

Images



All lanes : Rabbit IgG Isotype Control Rabbit mAb at 1/1k dilution
Lane 1 : Hela whole cell lysates
Lane 2 : 3T3 whole cell lysates
Lane 3 : C6 whole cell lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Exposure time: 2 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.