

Rat IL-1 β ELISA Kit

Catalog No: #EK0532

Package Size: #EK0532 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat IL-1 β ELISA Kit
Applications	ELISA
Species Reactivity	Rat

Application Details

Detect Range:62.5-4000 pg/mL

Sensitivity:15 pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Published Papers

Chenxing Wang;Siguang Ouyang;Xingjia Zhu;Yi Jiang;Zhichao Lu;Peipei Gong et al., Myricetin suppresses traumatic brain injury-induced inflammatory response via EGFR/AKT/STAT pathway, , (2023)

[PMID:38123650](https://pubmed.ncbi.nlm.nih.gov/38123650/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.