

Mouse (C3a) ELISA Kit

Catalog No: #EK18292

Package Size: #EK18292-1 48T #EK18292-2 96T



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse (C3a) ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:3ng/ml - 120ng/ml
Sample Type: Serum, plasma, Urine,tissue homogenates and other biological fluids.
Sample Volume:50ul
Assay Time: 3-5h
Detection wavelength: 618 nm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.