

Rat C3/Complement C3 ELISA Kit

Catalog No: #EK19388

Package Size: #EK19388 1x96T



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat C3/Complement C3 ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Specificity	Natural and recombinant Rat Complement C3
Target Name	Complement C3
Accession No.	P01026

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.072ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:C3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.