

Human C4BPB/C4b-binding protein beta chain ELISA Kit

Catalog No: #EK19903

Package Size: #EK19903 1x96T



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Human C4BPB/C4b-binding protein beta chain ELISA Kit
Applications	ELISA
Species Reactivity	Human
Specificity	Natural and recombinant Human C4b-binding protein beta chain
Target Name	C4b-binding protein beta chain
Accession No.	P20851

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.18ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:C4BPB

Published Papers

el at., Plasma proteome profiling reveals biomarkers of chemotherapy resistance in patients with advanced colorectal cancer, , (2024)
PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.