

Mouse Esm1/ Endothelial cell-specific molecule 1 ELISA Kit

Catalog No: #EK21129

Package Size: #EK21129-1 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Esm1/ Endothelial cell-specific molecule 1 ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:15.6-1000pg/mL

Sensitivity:8.02pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Detection Method:sandwich

Published Papers

Ni Kang;Yong Zhang;Fei Xue;Jinyu Duan;Fan Chen;Yu Cai;Qingxian Luan et al., Periodontitis induced by Porphyromonas gingivalis drives impaired glucose metabolism in mice, , (2022)

[PMID:36299624](https://pubmed.ncbi.nlm.nih.gov/36299624/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.