

EDIM ELISA Kit

Catalog No: #EK22008

Package Size: #EK22008 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EDIM ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range: Calculate Critical (CUT OFF): Critical = the average of Negative control well + 0.15. Negative Result: sample OD < Calculate Critical (CUT OFF) is Negative. Positive Result: sample OD > Calculate Critical (CUT OFF) is Positive.

Sample Type: Serum, plasma, Urine, tissue homogenates and other biological fluids.

Sample Volume: 50ul

Assay Time: 3-5h

Detection wavelength: 450 nm

Published Papers

Cen Li; Chen Xueyang; Lin Yiming; Liu Peihao; Mao Xinli; Ning Longgui; Shen Zhe; Xu Chengfu; Yu Chaohui; Yu Mengli; Zeng Hang; Zhang Hong; Zhang Jie; Zhang Yuwei et al., Vitamin D receptor targets hepatocyte nuclear factor 4 γ and mediates protective effects of vitamin D in nonalcoholic fatty liver disease, , (2020)

[PMID:32051143](#)

Zeng Zhang; Peng Wang; Tianyi Lei; Jianwei Guo; Yi Jiang; Yanhui Li; Jianxiong Zheng; Shunbing Wang; Haimuzi Xu; Guilin Jian; Quanbo Zhang; Yufeng Qing et al., The role and impact of the IL-6 mediated JAK2-STAT1/3 signaling pathway in the pathogenesis of gout., , (2025)

[PMID:40170729](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.