

APC anti-human CD23

Catalog No: #FC0029



Package Size: #FC0029-1 25Tests #FC0029-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description


|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Product Name       | APC anti-human CD23                                               |
| Host Species       | mouse                                                             |
| Clonality          | Monoclonal                                                        |
| Clone No.          | s-EBVCS-5                                                         |
| Isotype            | mouse IgG1,k                                                      |
| Applications       | FC                                                                |
| Species Reactivity | human                                                             |
| Conjugates         | APC                                                               |
| Target Name        | CD23                                                              |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with APC anti-human CD23 (colored histogram) or APC mouse IgG1 κ (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.