

APC anti-human/pig CD62E/CD62P

Catalog No: #FC0058



Package Size: #FC0058-1 25Tests #FC0058-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	APC anti-human/pig CD62E/CD62P
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-1.2B6
Isotype	mouse IgG1
Applications	FC
Species Reactivity	human/pig
Conjugates	APC
Target Name	CD62E/CD62P
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human Activated platelets stained with APC anti-human CD62E/CD62P (colored histogram) or PE-APC mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.