

Anti-L-Lactyl Lysine Rabbit mAb

Catalog No: #HL001

Package Size: #HL001-1 100ul



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-L-Lactyl Lysine Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	LR9001
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB, IHC-P, IP, ICC/IF, FCM, ChIP
Species Reactivity	Species independent
Target Name	L-Lactyl Lysine, Kla, L-la
Formulation	PBS pH 7.4
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC-P 1:200~1:1000 ICC/IF 1:50~1:100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.