

Anti-L-Lactyl Lysine Rabbit mAb

Catalog No: #HL001

Package Size: #HL001-1 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-L-Lactyl Lysine Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	LR9001
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB;IHC-P;IP;ICC/IF;FCM;ChIP
Species Reactivity	Species independent
Target Name	L-Lactyl Lysine, Kla, L-la
Formulation	PBS pH 7.4
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:1000 IHC-P 1:200~1:1000 ICC/IF 1:50~1:100

Images



All lanes : L-Lactyl Lysine Rabbit mAb at 1/1k dilution Lane 1 : Hela whole cell lysates Lane 2 : Hela treated with 100mM Sodium lactate for 24hr whole cell lysates Lysates/proteins at 20 μ g per lane. Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution Exposure time: 10 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.