

Histone H3 (L-Lactyl-Lys18) Rabbit mAb

Catalog No: #HL005

Package Size: #HL005-1 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histone H3 (L-Lactyl-Lys18) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	>95%. Protein A/G purified from cell culture supernatant.
Applications	WB;IHC;IF;ChIP
Species Reactivity	Human;Mouse;Rat
Conjugates	Unconjugated
Target Name	Histone H3
Other Names	H3/a, H3/c, H3/d, H3/f, H3/h
Accession No.	P68431
Calculated MW	15 kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20Λ C

Application Details

WB 1:500-1:2000; IHC 1:100-1:300; IF 1:100-1:300; ChIP 1:20-200;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.